

Stock up for spring! Athletic Tapes & Accessories



M-Tape® Zinc Oxide Trainers Tape

100% cotton backcloth, zinc oxide trainers tape. White, bleached M-Tape® features high tensile strength with excellent tackiness and conformability.

#130105 1 1/2" x 15 yd,
 32 rolls/case.....\$39.95

Mueller® Tear-Light Tape™

Lightweight elastic adhesive tape. Superior adhesive with excellent conformability. Easy to tear. Easy unwind. Ideal for securing lateral knee braces, as a compression wrap or fixation bandage for padding and cold packs.

#130626 1 1/2" x 7 1/2 yd, 32 rolls/case.....\$45.95
 #130632 2" x 7 1/2 yd, 24 rolls/case.....\$45.95



Sensi-Wrap (Blood Spill Tape)

- Cohesive stretch tape
- Adheres only to itself, not to patients' skin or hair.
- No damage to clothing or equipment
- Easy tear
- No scissors required
- Stops bleeding
- Apply directly to the wound area

#3172 2" x 5 yd,
 Individual.....\$1.55
 #3173 3" x 5 yd,
 Individual.....\$2.15



M-Wrap® (Underwrap) Natural

Pre-taping foam underwrap helps protect skin from tape chafing. Can be used to hold pads and socks in place and as a protective wrap inside hiking boots, ski boots, and other athletic footwear.

#130702 2 3/4" x 30 yd,
 48 rolls/case.....\$31.95



Tuffner™ Pre-Tape Spray

Spray adherent helps secure tapes and wraps.

#200902N 10 oz Aerosol.....\$6.95
 #200901N 4 oz Aerosol.....\$3.95

Tape Cutters and Shears



Double Ring Bandage Scissors

#S-5605 7 1/2".....\$6.95



Pro Trainers Scissors

#46360 8".....\$9.95



M-Cutter™

Athletic tape cutter cuts fast, safely. Features snap-in stainless steel blade cartridge. Replacement blades available.

#100201 cutter.....\$9.25
 #100202 replacement blade.....\$3.15

Stock up for spring! Bandages & Dressings



Woven Fabric Bands

Non-stick pad is highly absorbent and helps to cushion wound. Flexible adhesive fabric stretches and conforms to body contours for effective and comfortable protection. Latex free

- #016459 1" x 3", 100/Box.....\$3.65
- #013614 2" x 4 1/2", Extra Large, 50/Box.....\$5.95
- #016470 3" x 3", 4-Wing, 50/Box.....\$8.95
- #011350 1 3/4" x 2", Fingertip, 100/Box.....\$5.95
- #011325 1 3/4" x 3", Large Fingertip, 100/Box.....\$6.95
- #016425 Knuckle, 100/Box.....\$5.95



Non-Adhering Sterile Pads

The Non-Adherent Pads are designed to protect the wound while promoting the healing process. The non-stick outer facing will allow the dressing to absorb but not adhere to the wound. The easy removal of the pad will minimize the trauma to the wound site

- #3423 2" x 3", 100/Box.....\$6.95
- #3434 3" x 4", 100/Box.....\$8.95



Stretch Gauze Non-Sterile

This knitted stretch gauze bandage provides light compression to the site while securing holding the dressing in place. The rayon/poly blended product is specifically designed not to slip while in use. There are three packaging options available to meet all needs and applications. "Clean" wrap is a non-sterile product individually packaged.

- #3103 3" x 5 yd, 12/Bag.....\$4.45

Other sizes available at www.CFAMedical.com.

Click [Bandages & Dressings](#), then [Surgical Dressings](#) (pg. 2)

Foam-Felt-Podiatry



Podiatry Kit

This convenient Podiatry Kit contains:

- 8 Corn Pads (C2, 1/8" Felt)
- 8 Corn Pads (C2, 1/4" Felt)
- 6 Corn Pads (C38, 1/8" Felt)
- 6 Corn Pads (C38, 1/4" Felt)
- 8 Blister Pads (B39, 1/8" Felt)
- 8 Blister Pads (B39, 1/4" Felt)
- 6 Blister Pads (B41, 1/8" Felt)
- 6 Blister Pads (B41, 1/4" Felt)
- 4 Arch Pads (LA25, 1/8" Felt)
- 4 Arch Pads (LA25, 1/4" Felt)
- 4 Heel Lifts (H22, 1/8" Felt)
- 4 Heel Lifts (H22, 1/4" Felt)
- 4 Heel Spurs (H21, 1/8" Felt)
- 4 Heel Spurs (H21, 1/4" Felt)
- 4 Metatarsal Pads (M43, 1/8" Felt)
- 4 Metatarsal Pads (M43, 1/4" Felt)
- 4 Metatarsal Pads (M10, 1/8" Felt)
- 4 Metatarsal Pads (M10, 1/4" Felt)
- 6 Felt Heel and Lace Pads (2 1/2" x 4")
- 2 Foam Squares (4" x 4" x 1/8")
- 2 Foam Squares (4" x 4" x 1/4")
- 4 Moleskin Squares (4" x 4")

- #900517 Podiatry Kit.....\$39.95

Stock up for spring!
Spenco®



Spenco® Blister Kit

A compact 5 1/2" x 4 1/4" resealable pouch for convenient protection, prevention and relief from friction, pressure and blisters. Kit contains 6 – 3" x 1 1/2" and 1 - 3" x 5" Adhesive Knit, 6 – 1" squares 2ND SKIN® Dressing and 1 – 3" x 5" sheet of Pressure Pad ovals. Easy to carry in purse, backpack, tool box, car glove compartment or keep at home.

#49-106.....\$7.45



Spenco® Second Skin Squares

This revolutionary hydrogel dressing helps protect against blisters, pressure and friction. Helps bring cooling relief to stings, bites, poison ivy and blisters. Comes in three packages for a variety of applications. Dispenser Jar of 200 1-inch Squares.

#10-637 1" x 1", 200/Jar.....\$29.95



Spenco® Second Skin Circles

This revolutionary hydrogel dressing helps protect against blisters, pressure and friction. Helps bring cooling relief to stings, bites, poison ivy and blisters. Comes in three packages for a variety of applications. Dispenser Jar of 48 3-inch Circles.

#10-603 3", 48/Jar.....\$41.95



More Spenco® products available at www.CFAMedical.com.

Click [Foam-Felt-Podiatry / Spenco®](#), and go to page 5.

Blood Borne, Disinfectants & Gloves



Blood Clotting Spray

This revolutionary hydrogel dressing helps protect against blisters, pressure and friction. Helps bring cooling relief to stings, bites, poison ivy and blisters. Comes in three packages for a variety of applications. Dispenser Jar of 48 3-inch Circles.

#280540 3 oz.....\$5.95



Blood Stain Remover

Cleans and removes organic stains and odors. Its enzymes and surfactants rapidly convert water insoluble stains into water soluble stains. Biozyme is effective on blood, saliva, vomit, urine, fecal, wine, coffee and grass stains. (This product is not a disinfectant) 16 oz size.

#59089 16 oz.....\$8.95



Nasal Plugs

Soft, non-linting 100% virgin cotton plugs used to control common nosebleeds occurring in the first inch of the nose. 1 1/2" long

#191221 1 package of 300.....\$11.95

Stock up for spring!

Blood Borne, Disinfectants & Gloves (cont.)



Antimicrobial Hand Gel Wipe

A healthcare personnel handwash with 99.99% effective reduction against gram positive and gram negative bacteria, fungi, and yeast. Sani-Dex ALC physically removes bacteria and soil from hands – a function rub-in alcohol handwashes can not accomplish.

#D20300 100/Box.....\$5.95



Antiseptic Hand Gel

Extra thick AprilFresh® Antiseptic Gel quickly cleans and sanitizes hands and skin. It is ideal for situations where water is not readily available or extra protection is needed. A special skin conditioner prevents skin from drying out. 70% Ethyl Alcohol.

#11227 8 oz bottle
 w/ flip-top cap.....\$3.95



Sani-Cloth® Plus

Pre-moistened cloths contain powerful disinfectant to kill germs on hard surfaces and equipment. Convenient to use – economical, disposable. Tuberculocidal, bactericidal and virucidal including HIV. EPA registered.

#Q89072
 160 cloths/canister.....\$7.95



Biozide™ Disinfectant/Deodorant

Biozide™ is a ready-to-use, hospital grade, hard surface disinfectant. Our Phenolic formula is a fast acting, all-purpose, broad spectrum disinfectant/deodorant for environmental surfaces. Choosing Biozide™ will help you comply with the OSHA Bloodborne Pathogens Standard. Biozide™ is tuberculocidal, and virucidal. For use on: glass, porcelain, ceramic, metal, vinyl, polypropylene, polyester, rubber, bakelite, or any surfaces not damaged by alcohol. For use in: industrial, institutional, commercial, residential and medical facilities, equipment and vehicles. EPA Registered.

Disinfectant/Deodorant:

- Ready to use hospital grade disinfectant
- Kills HIV, TB and Herpes 1 & 2
- Virucidal, Tuberculocidal, Fungicidal, Bacterial
- Controls mold/mildew
- Neutralizes odors and freshens in one step
- Helps prevent cross-contamination
- EPA- registered

#34105 16 oz.....\$6.95



Mueller® WHIZZER® Cleaner & Disinfectant

KILLS STAPHYLOCOCCUS AU-REUS METHICILLIN-RESISTANT (MRSA)

Kills HIV-1, Herpes Simplex 1&2 & MRSA on hard, nonporous surfaces.

Versatile

Whizzer® is the perfect choice for homes, schools, locker rooms, camps, fitness centers, salons, daycare centers and hospitals. Anywhere people gather.

Use Whizzer® to kill bacteria on floors, walls, training tables, benches, door knobs, shower and bath areas, toilets, counters, athletic equipment, whirlpools and more! Whizzer is a concentrated cleaner that is effective against a broad spectrum of bacteria.

Powerful

Whizzer® kills HIV-1 (AIDS Virus), Hepatitis B & C, Herpes Simplex Virus Types 1 and 2, Norwalk Virus, SARS and other viruses. Kills strep, staph (including MRSA) and fungus on contact. Whizzer® helps prevent the spread of athlete's foot. It's a cleaner, fungicide, mildewstat, virucide and deodorizer all in one! EPA Registered.

Economical

One ounce of concentrate mixed with water makes a gallon of full-strength disinfectant. You get 128 gallons of disinfectant in every gallon of Whizzer®!

Two Free Spray Bottles with every Gallon!

EPA Registered

#230201 1 gallon/case (w/ 2 free sprayers).....\$34.95

Whizzer® available online at www.CFAMedical.com. Search "whizzer" in the search blank.

Stock up for spring!

Blood Borne, Disinfectants & Gloves (cont.)



Nitrile Exam Gloves (Non-Latex)

- The ideal solution for individuals sensitive to natural rubber latex and donning powder
- Contains no allergy causing natural rubber proteins
- Extraordinary strength and puncture resistance while maintaining tactile sensitivity

#A700-2511 Small 100/Box.....\$9.95
 #A700-2512 Medium 100/Box.....\$9.95
 #A700-2513 Large 100/Box.....\$9.95
 #A700-2514 XL 100/Box.....\$9.95



Vinyl Exam Gloves (Powder Free Non-Sterile)

- For unmatched quality and sensitivity, Dynarex vinyl exam gloves provide superior comfort, strength and stretch
- Pre-powdered for easy donning
- Ambidextrous

#A500-2611 Small 100/Box.....\$6.95
 #A500-2612 Medium 100/Box.....\$6.95
 #A500-2613 Large 100/Box.....\$6.95
 #A500-2614 XL 100/Box.....\$6.95

Braces & Supports



Mueller® Soft Ankle Brace with Straps

This lightweight, strapped ankle brace provides tape-like support without the hassle or expense of taping. The attached stirrup straps and figure 8 straps are less bulky than competitive braces, while offering comparable strength and support. The Soft Ankle Brace material and design conforms to either the left or right foot, providing comfort during competition.

- Secure stirrup straps provide support without being difficult to apply
- Imitate effects of taping the ankle with non-stretch stabilizing straps
- Durable 1000 denier polyester material provides strength and eliminates the need for metal eyelets
- Fits left or right foot

#213 Black (please specify size: XS, S, M, L, XL)\$19.95

Mueller Ankle Braces & Supports Sizing: (According to shoe size)

	XXS	XS	SM	MD	LG	XL
	<small>(Soccer only)</small>					
U.S. Men	4 - 5	6 - 7	7 - 9	9 - 11	11 - 13	13 - 15
U.S. Women	5 - 6	7 - 8	8 - 10	10 - 12	12 - 14	14 - 16

Mueller® High Performance ATF® Ankle Brace

Custom left or right fitted brace is designed to conform to the arch of the foot for precise fit and no bulk. Vertical stretch elastic and horizontal stretch elastic over the Achilles tendon allow complete plantar and dorsal flexion. Brace allows for unrestricted foot movement for jumping, speed, and agility. Specify left or right.

#M-201 Left (please specify size: XS, SM, MD, LG, XL)\$17.95
 #M-206 Right (please specify size: XS, SM, MD, LG, XL)\$17.95



AVAILABLE IN BLACK ONLY

Stock up for spring!

Braces & Supports (cont.)



Mueller® Jumper's Knee Strap

Mueller's Knee Strap features a tubular insert which provides focused mild pressure on the tendon below the knee cap to help guide the patella. Ideal relief for patellar tendonitis caused by frequent running and jumping, irritated kneecap and Osgood-Schlatter's Disease. Strap is fully adjustable and comfortable for all-day wear.

- #M-991 Red\$5.95
- #M-992 Black\$5.95
- #M-993 Blue\$5.95
- #M-994 Gold\$5.95

Mueller® Hinged Knee Brace

Patented triaxial hinge provide maximum support and protection and help prevent hyper-extension of the knee. Patella buttress stabilizes and protects the kneecaps. Opposing adjustable front and back straps hold brace in place. Contoured construction for comfortable fit minimizes bulk behind knee joint. Nylon lining for comfort. Fully trimmed edges for durability.

- #2333 Black (specify SM, M, LG, XL, XXL)\$34.95



Mueller® PRO-LEVEL™ Hinged Knee Brace Deluxe

Fits up to 28" thighs! Patented triaxial hinge produces near-normal joint motion- allow full mobility and provide maximum medial lateral support and protection in all sports or activities. Brace includes lock out option for complete immobilization and adjustable extension stops for limited range of motion if used as a rehab brace.

Features wraparound design for easy on/off, adjustable straps above and below the knees and six sizes SM-XXXL for a custom fit. The removable patella can be modified as needed and the open back helps eliminate bunching behind the knee. Extra-long (16") for maximum stability and protection.

- #B5333 Black (specify SM, M, LG, XL, XXL, XXXL)\$44.95

SIZING: (Measure around center of knee)		
SM	12 – 14"	30 – 35 cm
MD	14 – 16"	35 – 40 cm
LG	16 – 18"	40 – 45 cm
XL	18 – 20"	45 – 50 cm
XXL	20 – 22"	50 – 55 cm



Universal Shoulder Support



Universal Shoulder Support is an adjustable shoulder support with Velcro® closures. Outstanding features include: 4-way stretch neoprene for therapeutic heat retention to promote healing; Velcro® closures for adjustment at the arm and chest; and pockets over the AC joint to hold an ice pack, gel pack or extra padding. Also included is a

removable non-stretch Velcro strap to limit arm movement and restrict extension for chronic shoulder separation.

Fits either shoulder. Ideal for relief of sprains, strains, bursitis and tendinitis, and to control arm movement in the treatment of dislocations. *Measure around chest to determine size (see below).*

- #462 (specify M – up to 38", L – 38"-58", XL – 60"-68")\$49.95

Cold Shoulder Wrap



Cold therapy is one of the best ways to combat pain and swelling due to over exertion and heavy stress on muscles and joints. The Cramer Cold Shoulder Wrap provides complete cold therapy coverage for the shoulder, rotator cuff, upper arm and elbow to athletes who place a heavy workload on that area.

The Cold Shoulder Wrap is designed to be portable; one size fits most, right/left solution to icing athletes in and easy to use package. The sealed ice cells prevent leaking and the anti-microbial finish helps to prevent odor.

- #279826.....\$56.95

Hydration Products & Instruments

Lightning Detectors & More



SkyScan® Deluxe

Detects storms within 40 miles
 Easily tell if the storm is moving towards, away, or parallel to your location
 Identifies severe thunderstorms as well as storms that can produce dangerous winds, heavy rains or tornadoes
 Starts detecting thunderstorms up to 70 miles away, registers activity at 40 miles
 Used by little league teams, golf courses, and recreational districts
 Hand held or table mountable with optional wall mount accessory available
 Low battery indicator with battery save feature

Specifications:

- Ranges: 0–3 miles; 3–8 miles; 8–20 miles; and 20–40 miles
- Case: rugged, weather-resistant plastic case
- Power Source: 2 “9V” batteries, or optional 120V AC adapter
- Dimensions: 6.5" x 3.1" x 1.5"
- Weight: 10.8 oz (w/o batteries)

#LD50 SkyScan® Deluxe.....\$209.95



Compact Lightning Detector

Senses lightning within 6 miles, 6–12 miles, 12–24 miles and 20–40 miles. Power source: 1 “AAA” battery
 #LD9\$91.85



Digital Weather Radio

Mannix Digital NOAA Weather Alert Radio with audible storm alerts!

Features & Applications:
 Single line LCD displays N.O.A.A. alert, S.A.M.E. technology for geographic selection, Audible & visual alerts, Bright blue backlight, 12/24 hour clock with alarm, Day and date display, Digital PLL tuning, Digital volume adjustment, LCD Display.

Features:

- Single line LCD displays N.O.A.A. alert;
- S.A.M.E. technology allows for geographic selection
- Audible & visual alerts
- Bright blue backlight
- 12/24 hour clock with alarm
- Day and date display
- Digital PLL tuning
- Digital volume adjustment
- LCD Display

Supplied with batteries.

#WR102\$76.95

Facts About Lightning

Did You Know?....

Average lightning stroke is 6 miles long

Average thunderstorm is 6–10 miles wide

Average thunderstorm travels at a rate of 25 miles per hour

On average, thunder can be heard over a distance of 3–4 miles depending on humidity, terrain and other factors

Approximately 100,000 thunderstorms occur in the United States each year. Approximately 10% of all thunderstorms are severe enough to produce high winds, flash floods and tornadoes

Hydration Products & Instruments (cont.)

Heat & Humidity



MANNIX Digital Sling Psychrometer

Microprocessor technology moves the Mannix Digital Psychrometer far beyond old sling types. The Mannix Digital Psychrometer is a precision temperature & humidity meter with features only available on instruments five times the price. This unit is extremely accurate and gives you a super fast response.

A precision temperature and humidity instrument, it measures wet bulb and dry bulb temperature.

#SAM990DW\$99.95



MANNIX Digital Pocket Heat Index Monitor

Mannix Digital Pocket Heat Index Monitor displays heat index, temperature, humidity, & dew point! Ideal for monitoring factory & industrial conditions. Danger icon with audible beep, Min/max memory storage, Large digital display, F/C readings, Data hold, Reset feature, Auto power off.

#SAM800HI\$109.95

Hydration Testing

Portable Lighted Clinical Refractometer

Tri-scale portable Clinical Refractometer in a blue, foam-lined hard case. The ZGRC series of portable refractometers is designed to measure Urine Specific Gravity, Blood Plasma Serum Protein and Refractive Index, but includes the benefit of an internal light. The extra light can be turned on when the ambient light is so-so, and when light is good, you can turn off the light and use the unit!

Included are: a pipette for dropping test fluids onto the prism (never dip the refractometer into the test fluid!), a mini screwdriver for scale calibration, (3) 1.5 volt button batteries, and a Users Guide.


#ZGRC-200ATC\$129.95



Hydration Test Strips (Compare to Multistix® 10 SG Strips)

Tests for hydration. Tests urine for: Glucose. Bilirubin. Ketone. Specific gravity. Blood. pH. Protein. Urobilinogen. Nitrite. Leucocytes.

#H5701 (100/bottle).....\$49.95



Hydration & Cooling

Bottles Teampak™ (32)

Heavyweight, long-lasting sport drink bottles. Special plastic formula allows even pressure to dispense all liquid in the bottle. Wide mouth opening to easily add ice.

Economical, unassembled. Wide mouth quart bottles. Push/pull lid.

#020602 Quart (32/case) Color: Natural (clear).....\$64.95



Plastic Drinking Cups (10 oz)

#0050 10 oz. (1000/case).....\$33.95



Paper Products

Disposable Towels

3 ply
 #7227 13.5" x 18"
 (500/case).....\$21.95



Diagnostics



Economy Stethoscope Dual Head
 #900D.....\$6.95

Advantage™ 6013 Automatic BP Monitor

Smart Logic™ Technology for the ultimate in simplicity.

The 6013 features:

- Measures systolic, diastolic pressure, and pulse
- Smart Logic™ Technology determines ideal inflation level
- Automatic inflation and deflation control
- Self adjusting latex-free D-Bar cuff
- Recalls last 14 readings with date and time stamp
- Auto off

#6013-A Adult.....\$84.95
 #6013-L Large Adult.....\$94.95



Advantage™ 6016 Digital Wrist BP Monitor

The 6016 features:

- Measures systolic, diastolic pressure, and pulse
- New 8 bit microprocessor and advanced algorithm for more accurate readings
- Smart Logic™ Technology determines ideal inflation level
- Average mode automatically averages 3 readings for more meaningful results
- Latex-free wrist cuff fits 5" to 7" wrists
- Advanced 2 zone memory with date and time stamp recalls last 30 readings for two different people
- Auto off
- Compact storage case
- Optional thermal printer (sold separately)

#6016.....\$94.95



BP Cuff Economy Complete

Includes:

- Exclusive nylon cuff
- Latex or latex-free inflation bladder with two 55cm tubes
- Bulb and valve with filter screen protection
- For use with most aneroid or mercurial blood pressure instruments
- Individually boxed
- One Year Warranty
- Luer Connector included



#800-A Adult.....\$17.95
 #800-C Child.....\$18.95
 #800-X Large Adult.....\$26.95
 #800-T Thigh.....\$32.95

Detecto® Certifier Scale

A must for schools and universities with wrestling teams!

Detecto's 442 Certifier special purpose eye-level beam scale is "Legal for Trade." Many sporting events, such as wrestling as well as weight loss and fitness programs, require a scale which can be certified by weights and measures officials. The 442 Certifier fills this need.

This scale is designed and manufactured to maintain accuracy necessary to meet the requirements of Handbook No. 44, issued by the National Institute of Standards and Technology.

442 Specifications:
 Capacity: 350 lb x 4 oz or 160 kg x 100 g
 Platform Size: 13 3/4" x 18 1/2" / 349mm x 470mm
 Shipping Weight: 68 lb / 31 kg
 Color: White with gray mat
 Features: Dual wheels and non-slip mat

#442.....\$565.00





Diagnostics (cont.)

Detecto® Eye-Level Physician Scale

Detecto's model 439 features a height rod in inches and centimeters. The platform cover is easily removed for cleaning. The die-cast beam can be easily read from either side.

Detecto continues to set the pace in the clinical scale industry, backed by over a century of leadership. Every Detecto physician scale features a time-honored reputation for rigid construction and precision engineering. This assures an unbeatably rugged product that guarantees accuracy and durability. All Detecto scales are tested with weights traceable to the U.S. Bureau of Standards.

439 Specifications:
 Capacity: 400 lb x 4 oz
 Platform Size: 10 1/2" x 14 1/2" / 267 mm x 368 mm
 Shipping Weight: 38 lb / 17.2 kg
 Height Rod Measures: 30" – 78" and 76 cm – 198 cm
 Scale Height: 59" / 1.5 m

439 Features:

- Accuracy
- Heavy-duty understructure
- Die-cast beam, aluminum on black insert, easily read from either side
- Factory assembly
- Competitive pricing
- Detecto maintains the largest service network in the clinical scale industry

#D-439 without casters.....\$259.95



Skinfold Caliper

Lightweight aluminum construction. Measurement arms with "floating tips" under constant pressure assure accurate and reproducible results. Easy to read 60 mm scale. Booklet, tables and carrying case included.

#43063.....\$235.00



First Aid / Athletic Kits

The Roo™ Fanny Pack

Same features as the Roo™ Hand Warmers. Three zippered storage pockets. Main pocket contains 5 accessory loops for small items. Waterproof and wind proof. Measures 16"L x 7"H x 6"W. Belt expands to 58". Color black. Look for the C° ComfortTemp®.



#R-500.....\$29.95



Leather Athletic Trainer Holster

EMT Leather Holster Set Empty. Holds 7 instruments.

#5353.....\$29.95



First Aid Kit (Empty)

#112010 13" x 9 1/2" x 4 1/2".....\$14.95



Mueller® Sling Bag

Keep important items close at hand with this convenient shoulder sling bag that can easily be moved from back to front. The bag, with padded adjustable strap, can be worn over the shoulder or head. Hook and loop closures on pockets keep items from falling out. Exterior pockets allow easy access to key items: specifically made to fit rolls of tape, cell phone/pager and emergency shears/M-Cutter. Partitioned interior pockets and mesh pouches allow easy organization. Measures 12" x 4" x 14", central compartment capacity 396 cu. in., outer pocket capacity 42 cu. in.

#16017 12" x 4" x 14".....\$27.95



Cramer® Sling Pack

- Full zip compartments allow for easy unobstructed access to both main pockets.
- Exterior mesh pocket for quick access.
- Integrated cell phone holster also holds scissors for easy access.

#111985 15" x 10.5" x 5".....\$27.95



First Aid / Athletic Kits (cont.)



Trauma Bag ADC

Incredibly durable, extraordinarily versatile medical trauma bag for the pre-hospital care professional.

ADC Trauma Bag features:

- Designed for easy access and convenient storage.
- Constructed of 600 denier nylon for rugged durability
- Flap top with snap lock buckles permits rapid access
- Twelve elastic interior tool holders and three spacious padded compartments store needed equipment safely and securely
- Two zippered side pouches and dual front pleated pockets
- Includes an adjustable, removable webbed shoulder strap and contoured handle for easy carrying
- Zippered storage within the top flap
- Measures 17 1/2" x 10 3/4" x 7"
- Weighs 20 oz.
- One year warranty

#R-074.....\$54.95



Trauma Bag Neon Orange

Trauma Bag Neon Orange 600 Denier Cordura Nylon water proof, water resistant. Size: 17"x9"x7" with three padded compartments, quick snap buckles and adjustable shoulder strap.

#R-073.....\$39.95

Miscellaneous

Mouthguards

- Fast-fitting co-polymer resin
- Thickened posterior area helps absorb impact and oral trauma via an aerated shock absorbing "pillow"
- Impact pad reduces the risk of concussion
- Thin front wall construction does not interfere with speech or breathing
- Low profile helps reduce gagging
- Perfect fit gives better blunt-blow protection
- Co-polymer construction lends to comfort and durability
- Latex Free

Colors: Clear M-131007, White M-131032, Gold M-131013, Orange M-131035, Red M-131015, Maroon M-131033, Kelly Green M-131016, Forest Green M-131036, Royal Blue M-131017, Navy Blue M-131031, Purple M-131034, Black M-131018

SKU#s above (specify color)
 100/case.....\$34.95



Pharmaceuticals



Alcohol Wipes 200/box

#15-6818.....\$2.65



Antiseptic Wipes 100/box

#155201.....\$2.45



Ivy-Rid® Spray 3 oz.

Relief of poison ivy & oak as well as insect bites.

#850.....\$3.25



Bandage Spray 3 oz.

#151011.....\$3.65



Blood Clotting Spray 3 oz.

#280540.....\$5.95



Coolant Cold Spray 4 oz.

Non-medicated propellant to cool skin fast and help relieve swelling.

#030201.....\$3.45



Insect Sting Relief Wipes 50/Box

#157830.....\$5.95

Pharmaceuticals (cont.)



Triple Antibiotic Ointment 1 oz.

Generic antibiotic ointment prevents infection and aids healing of minor cuts, burns, and abrasions. Contains Polymyxin B Sulfate, Bacitracin Zinc and Neomycin Sulfate. 1 oz tube.

#906.....\$3.75



Triple Antibiotic Ointment 144/box

#907.....\$13.75



Hydrocortisone 1% 1 oz.

#203020.....\$2.45



Hydrocortisone 1% Unit Dose 144/box

Packaged: Unit dose, foil pack

#231409.....\$12.95



Save-A-Tooth® Emergency Tooth Preserving System

Protects and nourishes knocked-out teeth for replantation by a dentist. Uses a specially formulated pH balanced solution, as well as a patented unique removable basket and net container that protects tooth root cells crucial to the success of the replantation process.

#6302.....\$16.95



Ibuprofen Pain Relief Tablets Industrial Packets 250x2

Contains the same pain relief as two regular tablets of aspirin or acetaminophen. Sodium and sugar free. Contains 200mg. of Ibuprofen. Industrial packed 250 packets of 2.

#166185.....\$14.95



Non Aspirin 250x2

Contains Acetaminophen, the same ingredient as Tylenol®. Regular strength, 325 mg. tablets. Industrial packed 250 packets of 2.

#161581.....\$8.45



Electrolyte Heat Aid® 250x2

Heat Aid offers sodium free replacement of the body electrolyte lost through excessive perspiration, to help prevent heat fatigue and muscle cramps. It resupplies the body with the Potassium, Calcium, Magnesium, Phosphate and Chloride electrolytes taht are contained in perspiration. 250 packets of 2.

#21500.....\$12.95



Cramp Tablets 125x2

Low-cost relief of menstrual discomfort. Contains Acetaminophen and Pamabron. Sugar and salt free. 125 packets of 2 tablets per box. Package information presented in english and spanish.

#81048.....\$12.95



Diphen (Compare to Benadryl®) 200/box

Ingredients: Diphenhydramine HCL, 25MG. (Antihistamine). Packaged: 1 tablet per package. 200/box.

#18247.....\$10.95



Sinus Decongestant (Compare to Sudafed® PE) 500/box

Ingredients: Phenylephrine HCL 10 MG. Packaged: 1 tablet per packet. 500/box

#26360.....\$21.50



Coradin Forte 500/box

For cold, allergy and flu symptoms. Ingredients: Acetaminophen (Analgesic) 162.5 MG, Chlorpheniramine Maleate (Antihistamine) 2 MG, and Phenylephrine HCL (Decongestant) 5 MG. Packaged: 2 tablets per packet. 500/Box.

#26356.....\$16.95



Antacid (Compare to Tums®) 500/box

Ingredients: Calcium Carbonate 420 MG. Packaged: 2 Tablets per package. 500/box.

#171547.....\$15.95



Diotame (Compare to Pepto Bismol®) 100/box

Antacid/Diarrhea relief. Ingredients: Bismuth Subsalicylate 262 MG. Packaged: 2 tablets per package. 100/box.

#22033.....\$8.95

Rehab

More Rehab products available at www.CFAMedical.com. Click [Rehab](#), [Hydrocollators](#), [Muscle Stims](#), & [Whirlpools](#)

Intellect® Transport Electrotherapy

'On-the-Go' Electrotherapy

- 2 channels of electrical stimulation output
- Ten user-defined memory positions for User Protocols
- 110V or battery-powered option
- Independent intensity and parameter controls for each channel
- 4 standard waveforms on board: Interferential, Premod, Russian and High Volt
- Logical control system and a large easy to read graphical LCD display
- Unit uniquely designed for tabletop, wall mount, or mobile use.



#2783.....\$1645.00

Intellect® Transport Ultrasound

A New Dimension in Ultrasound

- Fully functional 1 and 3.3 MHz frequencies
- Ten user-defined memory positions for User Protocols
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Beam Nonuniformity Ratio (BNR) 5:1 maximum
- 110V or battery powered
- Unit uniquely designed for tabletop, wall mount, or mobile use
- Ergonomically designed sound head
- Applicators in variety of sizes (1 cm², 2 cm², 5 cm² and 10 cm²)
- All sound head applicators are interchangeable due to Electronic Signature™ residing on all transducers
- Head warming feature on all sound head applicators



#2782.....\$1195.00

Hydrocollator® E-1 Stationary Heating Unit

- High quality stainless steel
- Thermostatically controlled temperature
- No plumbing required



Product Weight : 15 lbs (7 kg)
 Product Length : 13 in (33 cm)
 Product Width : 8 in (20 cm)
 Product Height : 16 in (41 cm)

Includes 4 standard hot packs.

#2102.....\$315.00

LVES 1000

Low volt, low frequency electrical stimulator penetrates into deep tissue layers. Uses include: pain reduction, reduce inflammation, non-invasive acupuncture, enhance soft tissue healing, decrease edema, increase blood flow, fibromyalgia, lupus, reduce muscle tension causing headaches, etc. Great with heat and cold therapy. Small size provides excellent portability on the road and in the medical facility. Ten different treatment settings: 6 individual and 4 combination settings. Treatment time ranges 5 to 30 minutes. Intensity settings: 1 to 10. Uses 3 AAA batteries; included.



provides excellent portability on the road and in the medical facility. Ten different treatment settings: 6 individual and 4 combination settings. Treatment time ranges 5 to 30 minutes. Intensity settings: 1 to 10. Uses 3 AAA batteries; included.

#LVES 1000.....\$189.95

Splints

Visit our wide array of splinting products available at www.CFAMedical.com. Click [Splinting & Crutches](#).

Sam Splint

Made of flexible aluminum with foam padding.

#1121 36" x 4.25" roll.....\$15.95



Thermometry



ADTEMP I™ 412 Thermometer

The ideal mercury replacement thermometer.

Adtemp I™ Thermometer features:

- Super low cost digital thermometer
- Compact size – under 5" long, weighs just 14 grams
- Auto off function conserves battery life
- Range 90°F-109.9°F ±.2°F / 32°C – 43.9°C ±.1°C
- Non-replaceable 1.55v (LR41) type battery provides up to 2,500 measurements
- Integral case keeps instrument sanitary
- Beeps when ready

#412.....\$4.20



ADTEMP Digital Therm Sheaths 100/Box

Thermometer Sheaths feature:

- Non-sterile, disposable design helps eliminate cross contamination
- Cost effective & easy-to-apply
- Suitable for oral or rectal use
- Digital model designed for use with ADC and most comparably sized digital thermometers
- For use on all digital thermometers
- Made in the U.S.A.

#416-100 100/Box.....\$4.20

We've Only Scratched the Surface.

This Spring Sports Guide contains just a small selection of our huge inventory. For more great products, visit our website at

www.cfamedical.com

