

maximum

SUPPORT LEVEL

ANKLE

FITS LEFT OR RIGHT

THE ONE Ankle Brace

Anti-slip side straps help protect against inversion & eversion ankle sprains



Neoprene-free, Latex-free

Adjustable straps simulate figure 8 taping structure

Lasts much longer than tape, won't break down

For superior support this brace is *THE ONE*

Protected by the
aegis
MICROBE SHIELD®

Controls or eliminates odors, unsightly stains and product deterioration caused by the growth of bacteria and microorganisms.

El Uno Tobillo



MUELLER®
SPORT CARE®



maximum

SUPPORT LEVEL

**MUELLER[®]
SPORT CARE[®]**

FITS LEFT OR RIGHT

We equip athletes with the most advanced sports medicine products available.

STAY IN THE GAME™

Anti-slip Criss-Crossing

Side Straps help protect against inversion and eversion ankle sprains. Adjust them at any time to help to prevent brace slippage ensuring a tight, comfortable fit.



**THE ONE
Ankle Brace**



Top Elastic Finishing Straps lock down all of the straps and provide a clean finish



Side-Pulling Stirrup Straps imitate the effects of taping



Cushioned, breathable tongue keeps this brace extremely comfortable

If you need superior support, this ankle brace is *THE ONE*...

Very light and very low profile, *THE ONE* is totally comfortable to wear in shoes or cleats during competition. All the traditional benefits of taping, but with all the adjustable flexibility of a soft strap brace.

***THE ONE* is ideal for basketball, football, baseball, softball, volleyball, soccer, and other high-impact sports** to help support strained, sprained, twisted or swollen ankles.



Contact Your Authorized Mueller Team Dealer

www.MuellerSportsMed.com • 800-356-9522

FITS STANDARD SHOE SIZES:

| | XS | SM | MD | LG | XL | XXL | XXXL |
|-----------------------|---------|---------|----------|-----------|-----------|-----------|-----------|
| U.S. MEN | 6-7 | 7-9 | 9-11 | 11-13 | 13-15 | 15-17 | 17-19 |
| U.S. WOMEN | 7-8 | 8-10 | 10-12 | 12-14 | 14-16 | 16-18 | 18-20 |
| U.K. MEN | 5-6 | 6-8 | 8-10 | 10-12 | 12-14 | 14-16 | 16-18 |
| U.K. WOMEN | 5.5-6.5 | 6.5-8.5 | 8.5-10.5 | 10.5-12.5 | 12.5-14.5 | 14.5-16.5 | 16.5-18.5 |
| EUROPEAN MEN | 38-40 | 40-42 | 42-45 | 45-47 | 47-49 | 49-51 | 51-53 |
| EUROPEAN WOMEN | 36-39 | 39-41 | 41-43 | 44-46 | 46-48 | 48-50 | 50-52 |