

maximum

SUPPORT LEVEL

ANKLE

FITS LEFT OR RIGHT

THE ONE Ankle Brace

Anti-slip side straps help protect against inversion & eversion ankle sprains



Neoprene-free, Latex-free

Adjustable straps simulate figure 8 taping structure

Lasts much longer than tape, won't break down

For superior support this brace is *THE ONE*

Protected by the
aegis
MICROBE SHIELD®

Controls or eliminates odors, unsightly stains and product deterioration caused by the growth of bacteria and microorganisms.

El Uno Tobillo



MUELLER®
SPORT CARE®



maximum

SUPPORT LEVEL

**MUELLER[®]
SPORT CARE[®]**

FITS LEFT OR RIGHT

We equip athletes with the most advanced sports medicine products available.

STAY IN THE GAME™

Anti-slip Criss-Crossing

Side Straps help protect against inversion and eversion ankle sprains. Adjust them at any time to help to prevent brace slippage ensuring a tight, comfortable fit.

**THE ONE
Ankle Brace**



Top Elastic Finishing Straps lock down all of the straps and provide a clean finish



Side-Pulling Stirrup Straps imitate the effects of taping



Cushioned, breathable tongue keeps this brace extremely comfortable

If you need superior support, this ankle brace is *THE ONE*...

Very light and very low profile, *THE ONE* is totally comfortable to wear in shoes or cleats during competition. All the traditional benefits of taping, but with all the adjustable flexibility of a soft strap brace.

***THE ONE* is ideal for basketball, football, baseball, softball, volleyball, soccer, and other high-impact sports** to help support strained, sprained, twisted or swollen ankles.



Contact Your Authorized Mueller Team Dealer

www.MuellerSportsMed.com • **800-356-9522**

FITS STANDARD SHOE SIZES:

	XS	SM	MD	LG	XL	XXL	XXXL
U.S. MEN	6-7	7-9	9-11	11-13	13-15	15-17	17-19
U.S. WOMEN	7-8	8-10	10-12	12-14	14-16	16-18	18-20
U.K. MEN	5-6	6-8	8-10	10-12	12-14	14-16	16-18
U.K. WOMEN	5.5-6.5	6.5-8.5	8.5-10.5	10.5-12.5	12.5-14.5	14.5-16.5	16.5-18.5
EUROPEAN MEN	38-40	40-42	42-45	45-47	47-49	49-51	51-53
EUROPEAN WOMEN	36-39	39-41	41-43	44-46	46-48	48-50	50-52